

MK-3CH HEMATOLOGY CONTROL+CRP 2.5mL

ASSAY SHEET

LOT	MK3C248H
	2024-11-05
OPEN	14DAYS

MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.


1



2



3



		MEK-1303	RANGE
WBC	10 ³ /μL	20.35	± 2.20
RBC	10 ⁶ /μL	5.87	± 0.25
HGB	g/dL	19.32	± 0.70
HCT	%	62.9	± 4.0
MCV	fL	107.1	± 7.0
MCH	pg	32.9	± 2.8
MCHC	g/dL	30.7	± 3.0
PLT	10 ³ /μL	541.5	± 70.0
MPV	fL	8.6	± 3.0
RDW-CV	%	16.5	± 3.0
LY%	%	13.76	± 10.00
MO%	%	2.57	± 5.00
GR%	%	83.67	± 10.00
LY	10 ³ /μL	2.80	± 1.80
MO	10 ³ /μL	0.52	± 0.70
GR	10 ³ /μL	17.03	± 2.00
CRP	Precision	≤ 10.0% (52.90-98.30 mg/L)	

Manual Registration	
Lot	248
Level	High
Control	MK-3C



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN

CE 0123 IVD
0694-901767G

NIHON KOHDEN

MK-3CL HEMATOLOGY CONTROL+CRP 2.5mL

ASSAY SHEET

LOT	MK3C248L
	2024-11-05
OPEN	14DAYS

MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.


1



2



3



		MEK-1303	RANGE
WBC	10 ³ /μL	2.01	± 0.60
RBC	10 ⁶ /μL	2.27	± 0.15
HGB	g/dL	6.02	± 0.40
HCT	%	20.2	± 2.0
MCV	fL	88.8	± 9.0
MCH	pg	26.5	± 2.4
MCHC	g/dL	29.8	± 3.0
PLT	10 ³ /μL	73.3	± 30.0
MPV	fL	8.6	± 3.0
RDW-CV	%	15.7	± 3.0
LY%	%	66.31	± 10.00
MO%	%	5.85	± 5.00
GR%	%	27.84	± 10.00
LY	10 ³ /μL	1.33	± 0.50
MO	10 ³ /μL	0.12	± 0.20
GR	10 ³ /μL	0.56	± 0.50
CRP	Precision	≤ 10.0%	
		(1.30-7.30 mg/L)	

Manual Registration	
Lot	248
Level	Low
Control	MK-3C



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN

CE **IVD**
 0123
 0694-901765G

NIHON KOHDEN

MK-3CN HEMATOLOGY CONTROL+CRP 2.5mL

ASSAY SHEET

LOT	MK3C248N
	2024-11-05
OPEN	14DAYS

MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 ³ /μL	7.82	± 0.80
RBC	10 ⁶ /μL	4.79	± 0.20
HGB	g/dL	14.14	± 0.50
HCT	%	47.2	± 3.0
MCV	fL	98.5	± 7.0
MCH	pg	29.5	± 2.4
MCHC	g/dL	30.0	± 3.0
PLT	10 ³ /μL	266.6	± 50.0
MPV	fL	8.5	± 2.0
RDW-CV	%	16.1	± 3.0
LY%	%	30.78	± 10.00
MO%	%	3.40	± 5.00
GR%	%	65.82	± 10.00
LY	10 ³ /μL	2.41	± 1.00
MO	10 ³ /μL	0.27	± 0.40
GR	10 ³ /μL	5.14	± 1.00
CRP	Precision	≤ 10.0% (23.00-42.80 mg/L)	

Manual Registration	
Lot	248
Level	Normal
Control	MK-3C



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901766G