

**NIHON KOHDEN**

**MEK-CAL** HEMATOLOGY CALIBRATOR 2.0mL

**ASSAY SHEET**

LOT	PLUS24Y
	2024-12-05
OPEN	7DAYS

MEK-9100,9200  
2D Barcode Entry



MEK-1301,1302,1303,1305  
Barcode Entry



		MEK -6400	MEK -6500	MEK -7222	MEK -7300	MEK -8222	MEK -9100	MEK -1301
		-6410	-6510				-9200	-1302
		-6420	-6550					-1303
		-6450						-1305
WBC	10 <sup>3</sup> /μL	9.0	9.0	9.2	9.1	9.0	9.39	8.98
RBC	10 <sup>6</sup> /μL	4.64	4.64	4.59	4.63	4.59	4.58	4.62
HGB	g/dL	13.4	13.6	13.3	13.3	13.2	13.55	13.53
HCT	%	39.7	39.7	41.2	40.0	39.0	44.3	44.8
MCV	fL	85.6	85.6	89.7	86.5	84.9	96.8	96.9
MCH	pg	28.9	29.3	29.0	28.7	28.8	29.6	29.3
MCHC	g/dL	33.8	34.3	32.3	33.3	33.8	30.6	30.2
PLT	10 <sup>3</sup> /μL	240	240	235	238	237	250.4	258.3
MPV	fL	6.5	8.6	6.5	7.5	6.5	8.9	9.0
RDW-CV*	%	17.0	17.0	15.5	17.1	15.9	17.5	17.4

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	170
FL	87
SD	112
FS THR.	100

\* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN



0694-901780U