

**MK-3CN** HEMATOLOGY CONTROL+CRP 2.5mL

**ASSAY SHEET**

<b>LOT</b>	<b>MK3C24YN</b>
	<b>2025-02-05</b>
<b>OPEN</b>	<b>14DAYS</b>

**MEK-1303 Barcode Entry**

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 <sup>3</sup> /μL	7.46	± 0.80
RBC	10 <sup>6</sup> /μL	4.76	± 0.20
HGB	g/dL	13.78	± 0.50
HCT	%	45.6	± 3.0
MCV	fL	95.8	± 9.0
MCH	pg	28.9	± 2.4
MCHC	g/dL	30.2	± 3.0
PLT	10 <sup>3</sup> /μL	276.7	± 50.0
MPV	fL	8.5	± 2.0
RDW-CV	%	16.7	± 3.0
LY%	%	35.43	± 10.00
MO%	%	2.11	± 5.00
GR%	%	62.46	± 10.00
LY	10 <sup>3</sup> /μL	2.64	± 1.00
MO	10 <sup>3</sup> /μL	0.16	± 0.40
GR	10 <sup>3</sup> /μL	4.66	± 1.00
CRP	Precision	≤ 10.0% (20.90-38.90 mg/L)	

Manual Registration	
Lot	24Y
Level	Normal
Control	MK-3C



<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN



0694-901766H

**MK-3CH** HEMATOLOGY CONTROL+CRP 2.5mL

**ASSAY SHEET**

<b>LOT</b>	<b>MK3C24YH</b>
	<b>2025-02-05</b>
<b>OPEN</b>	<b>14DAYS</b>

**MEK-1303 Barcode Entry**

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 <sup>3</sup> /μL	19.80	± 2.20
RBC	10 <sup>6</sup> /μL	5.95	± 0.25
HGB	g/dL	19.17	± 0.70
HCT	%	62.7	± 4.0
MCV	fL	105.4	± 9.0
MCH	pg	32.2	± 2.8
MCHC	g/dL	30.6	± 3.0
PLT	10 <sup>3</sup> /μL	556.1	± 70.0
MPV	fL	8.5	± 3.0
RDW-CV	%	17.1	± 3.0
LY%	%	17.74	± 10.00
MO%	%	1.29	± 5.00
GR%	%	80.97	± 10.00
LY	10 <sup>3</sup> /μL	3.51	± 1.80
MO	10 <sup>3</sup> /μL	0.26	± 0.70
GR	10 <sup>3</sup> /μL	16.03	± 2.00
CRP	Precision	≤ 10.0% (50.20-93.20 mg/L)	

Manual Registration	
Lot	24Y
Level	High
Control	MK-3C



<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN



0694-901767H

# NIHON KOHDEN

**MK-3CL** HEMATOLOGY CONTROL+CRP 2.5mL

## ASSAY SHEET

<b>LOT</b>	<b>MK3C24YL</b>
	<b>2025-02-05</b>
<b>OPEN</b>	<b>14DAYS</b>

**MEK-1303 Barcode Entry**

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 <sup>3</sup> /μL	2.03	± 0.60
RBC	10 <sup>6</sup> /μL	2.40	± 0.15
HGB	g/dL	6.09	± 0.40
HCT	%	20.4	± 2.0
MCV	fL	85.0	± 9.0
MCH	pg	25.4	± 2.4
MCHC	g/dL	29.9	± 3.0
PLT	10 <sup>3</sup> /μL	69.8	± 30.0
MPV	fL	8.2	± 3.0
RDW-CV	%	15.6	± 3.0
LY%	%	64.29	± 10.00
MO%	%	3.70	± 5.00
GR%	%	32.01	± 10.00
LY	10 <sup>3</sup> /μL	1.31	± 0.50
MO	10 <sup>3</sup> /μL	0.08	± 0.20
GR	10 <sup>3</sup> /μL	0.64	± 0.50
CRP	Precision	≤ 10.0%	
		(1.70-9.90 mg/L)	

Manual Registration	
Lot	24Y
Level	Low
Control	MK-3C



<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**  
 1-31-4 NISHIOCHIAI, SHINJUKU-KU  
 TOKYO 161-8560 JAPAN



0694-901765H