


MEK-3DN HEMATOLOGY CONTROL 2.0mL

ASSAY SHEET

LOT	B252N
	2025-05-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order.
If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±0.8	8.0	8.0	7.88
RBC	10 ⁶ /μL	±0.20	4.90	4.90	4.93
HGB	g/dL	±0.5	14.3	14.5	14.34
HCT	%	±3.0	41.8	41.8	48.0
MCV	fL	±7.0	85.3	85.3	97.4
MCH	pg	±2.4	29.2	29.6	29.1
MCHC	g/dL	±3.0	34.2	34.7	29.9
PLT	10 ³ /μL	±50	252	252	269.1
MPV	fL	±2.0	5.3	8.0	8.2
RDW-CV *	%	±3.0	16.7	16.7	16.7
LY%	%	±10.0	31.8	31.8	28.74
MO%	%	±5.0	3.9	3.9	4.14
GR%	%	±10.0	64.3	64.3	67.12
LY	10 ³ /μL	±1.0	2.5	2.5	2.26
MO	10 ³ /μL	±0.4	0.3	0.3	0.33
GR	10 ³ /μL	±1.0	5.2	5.2	5.29

MEK-1305
QC ESR
7±5 mm

Manual Registration	
Lot	252
Level	Normal
Control	MEK-3D

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>



NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN

