

NIHON KOHDEN

MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS251
	2025-02-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



		MEK -6400	MEK -6500	MEK -7222	MEK -7300	MEK -8222	MEK -9100	MEK -1301
		-6410	-6510				-9200	-1302
		-6420	-6550					-1303
		-6450						-1305
WBC	10 ³ /μL	8.8	8.8	9.1	9.0	8.8	9.13	8.78
RBC	10 ⁶ /μL	4.56	4.56	4.55	4.59	4.55	4.54	4.58
HGB	g/dL	13.6	13.7	13.5	13.4	13.4	13.76	13.67
HCT	%	40.1	40.1	42.1	40.9	39.8	45.6	46.1
MCV	fL	87.9	87.9	92.5	89.1	87.5	100.5	100.7
MCH	pg	29.8	30.0	29.7	29.2	29.5	30.3	29.8
MCHC	g/dL	33.9	34.2	32.1	32.8	33.7	30.2	29.7
PLT	10 ³ /μL	224	224	225	228	227	240.7	248.3
MPV	fL	5.6	8.2	5.6	7.1	5.6	8.5	8.6
RDW-CV *	%	16.1	16.1	15.0	16.3	15.3	16.3	16.2

MEK	
	-7222
	-7300
	-8222
	-9100
	-9200
FS	168
FL	88
SD	110
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



MEK-CAL



PLUS251



2025-02-05



0694-901780W