

LOT	B255N
	2025-08-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order.

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±0.8	8.0	8.0	7.82
RBC	10 ⁶ /μL	±0.20	4.89	4.89	4.91
HGB	g/dL	±0.5	14.1	14.2	14.11
HCT	%	±3.0	41.6	41.6	47.5
MCV	fL	±7.0	85.0	85.0	96.7
MCH	pg	±2.4	28.8	29.0	28.7
MCHC	g/dL	±3.0	33.9	34.1	29.7
PLT	10 ³ /μL	±50	251	251	274.6
MPV	fL	±2.0	5.3	8.0	8.2
RDW-CV *	%	±3.0	16.7	16.7	17.1
LY%	%	±10.0	33.0	33.0	30.49
MO%	%	±5.0	3.4	3.4	2.55
GR%	%	±10.0	63.6	63.6	66.96
LY	10 ³ /μL	±1.0	2.6	2.6	2.38
MO	10 ³ /μL	±0.4	0.3	0.3	0.20
GR	10 ³ /μL	±1.0	5.1	5.1	5.24

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.

MEK-1305
QC ESR
8±5 mm

Manual Registration	
Lot	255
Level	Normal
Control	MEK-3D



<https://www.nihonkohden.com/spi/index.html>



NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



MEK-3DN



B255N



2025-08-05



0694-901685T

