

MEK-5DH

HEMATOLOGY CONTROL
for WBC 5 part differential

3.0mL

ASSAY SHEET

LOT	MEK5D255H
	2025-08-05
OPEN	14DAYS

MEK-9100,9200
2D Barcode Entry



		RANGE	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 ³ /μL	±2.2	25.0	24.9	24.4	25.63
NE%	%	±10.0	62.7	62.7	62.7	62.88
LY%	%	±10.0	20.9	20.9	20.9	23.11
MO%	%	±5.0	11.5	11.5	11.5	9.18
EO%	%	±5.0	3.5	3.5	3.5	3.68
BA%	%	±4.0	1.4	1.4	1.4	1.15
NE	10 ³ /μL	±2.8	15.6	15.6	15.3	16.13
LY	10 ³ /μL	±2.8	5.2	5.2	5.1	5.92
MO	10 ³ /μL	±1.0	2.9	2.9	2.8	2.35
EO	10 ³ /μL	±1.0	0.9	0.9	0.9	0.94
BA	10 ³ /μL	±0.8	0.4	0.3	0.3	0.29
RBC	10 ⁶ /μL	±0.25	5.27	5.24	5.23	5.27
HGB	g/dL	±0.7	16.5	16.2	16.3	16.28
HCT	%	±4.0	52.6	51.1	49.7	57.8
MCV	fL	±8.0	99.8	97.6	95.1	109.7
MCH	pg	±2.8	31.3	30.9	31.2	30.9
MCHC	g/dL	±3.0	31.4	31.7	32.8	28.2
RDW-CV *	%	±4.0	13.8	14.6	14.0	15.1
PLT	10 ³ /μL	±70	563	554	562	576.8
MPV	fL	±4.0	7.1	7.7	7.2	9.7

Manual Registration	
Lot	255
Level	High
Control	MEK5D

* RDW-CV is displayed as RDW on MEK-7222,8222.



<https://www.nihonkohden.com/spi/index.html>



NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901692S

Morelia (matriz) Tel: 443 690 2760
León (suc.) Tel: 477 510 2573
San Luis Potosí (suc.) Tel: 444 535 8886
www.multiservicioslab.com

