

LOT	B258L
	2025-11-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order.
If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±0.6	2.1	2.1	2.04
RBC	10 ⁶ /μL	±0.15	2.33	2.33	2.26
HGB	g/dL	±0.4	5.8	6.0	6.02
HCT	%	±2.0	18.2	18.2	20.2
MCV	fL	±9.0	78.3	78.3	89.6
MCH	pg	±2.4	24.9	25.8	26.6
MCHC	g/dL	±3.0	31.9	33.0	29.8
PLT	10 ³ /μL	±30	62	62	70.5
MPV	fL	±3.0	6.3	8.6	8.3
RDW-CV *	%	±3.0	17.4	17.4	15.5
LY%	%	±10.0	60.6	60.6	60.90
MO%	%	±5.0	6.9	6.9	3.09
GR%	%	±10.0	32.5	32.5	36.01
LY	10 ³ /μL	±0.5	1.3	1.3	1.24
MO	10 ³ /μL	±0.2	0.1	0.1	0.06
GR	10 ³ /μL	±0.5	0.7	0.7	0.74

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.

MEK-1305
QC ESR
50±17 mm

Manual Registration	
Lot	258
Level	Low
Control	MEK-3D



<https://www.nihonkohden.com/spi/index.html>



NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN

