

MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS25Y
	2025-12-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



MEK-CAL



<https://www.nihonkohden.com/spi/index.html>



PLUS25Y



2025-12-05



0694-901832F

		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	9.2	9.2	9.2	9.1	9.34	9.10
RBC	10 ⁶ /μL	4.53	4.53	4.51	4.46	4.44	4.45
HGB	g/dL	13.0	13.2	13.0	12.9	13.22	13.15
HCT	%	39.1	39.1	39.6	38.5	44.0	44.2
MCV	fL	86.3	86.3	87.7	86.3	99.0	99.4
MCH	pg	28.7	29.1	28.8	28.9	29.8	29.6
MCHC	g/dL	33.2	33.8	32.8	33.5	30.0	29.8
PLT	10 ³ /μL	218	218	218	218	225.1	231.8
MPV	fL	5.7	8.2	7.1	5.7	8.2	8.3
RDW-CV *	%	15.9	15.9	16.1	15.1	15.9	15.9

MEK	-7300 -8222 -9100 -9200
FS	171
FL	88
SD	109
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,8222.



NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN