

MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS263
	2026-04-05
OPEN	7DAYS

		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	8.8	8.8	8.9	8.7	9.05	8.83
RBC	10 ⁶ /μL	4.62	4.62	4.61	4.57	4.56	4.62
HGB	g/dL	13.5	13.7	13.5	13.4	13.76	13.73
HCT	%	40.7	40.7	41.0	39.9	45.7	46.6
MCV	fL	88.2	88.2	89.0	87.4	100.3	100.8
MCH	pg	29.2	29.7	29.3	29.3	30.2	29.7
MCHC	g/dL	33.2	33.7	32.9	33.6	30.1	29.5
PLT	10 ³ /μL	220	220	220	219	230.1	240.9
MPV	fL	6.0	8.4	7.2	6.0	8.5	8.5
RDW-CV *	%	16.0	16.0	16.0	15.0	15.8	15.8

MEK	
	-7300
	-8222
	-9100
	-9200
FS	163
FL	88
SD	99
FS THR.	100

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901832J