

MEK-5DN HEMATOLOGY CONTROL for WBC 5 part differential 3.0mL

ASSAY SHEET

LOT	MEK5D265N
	2026-08-05
OPEN	14DAYS

MEK-9100,9200
2D Barcode Entry



		RANGE	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 ³ /μL	±0.8	8.2	8.0	8.21
NE%	%	±15.0	63.4	63.4	64.00
LY%	%	±15.0	22.1	22.1	24.06
MO%	%	±5.0	10.8	10.8	7.98
EO%	%	±5.0	2.4	2.4	2.70
BA%	%	±4.0	1.3	1.3	1.26
IG% **	%	±4.0	-	-	0.54
NE	10 ³ /μL	±1.5	5.2	5.0	5.25
LY	10 ³ /μL	±1.5	1.8	1.8	1.98
MO	10 ³ /μL	±0.5	0.9	0.9	0.66
EO	10 ³ /μL	±0.5	0.2	0.2	0.22
BA	10 ³ /μL	±0.4	0.1	0.1	0.10
IG **	10 ³ /μL	±0.4	-	-	0.04
RBC	10 ⁶ /μL	±0.20	4.40	4.45	4.45
HGB	g/dL	±0.5	13.1	13.0	13.36
HCT	%	±3.0	39.7	39.1	45.0
MCV	fL	±6.0	90.3	87.8	101.1
MCH	pg	±2.8	29.8	29.2	30.0
MCHC	g/dL	±3.0	33.0	33.2	29.7
RDW-CV *	%	±4.0	15.2	14.4	14.8
PLT	10 ³ /μL	±50	266	245	254.7
MPV	fL	±2.0	5.5	4.2	7.0

Manual Registration	
Lot	265
Level	Normal
Control	MEK5D

* RDW-CV is displayed as RDW on MEK-8222.

** MEK-9200 only
Only software versions 02-01 and later can be used.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



MEK-5DN



265N



2026-08-05



0694-901691W