



Instruments:		Norma iVet-5																	
Parameter		CONTROL "Low"				CONTROL "Normal"				CONTROL "High"									
		LOT: 30561011		LOT: 30561012		LOT: 30561013		LOT: 30561013											
		Mean	±	Limit	Range		Mean	±	Limit	Range		Mean	±	Limit	Range				
					Variation					Variation					Variation				
WBC	10 ⁹ /L	3.14	±	0.50	2.64	-	3.64	8.08	±	1.00	7.07	-	9.07	15.91	±	2.50	13.41	-	18.41
RBC	10 ¹² /L	2.45	±	0.25	2.20	-	2.70	4.43	±	0.30	4.13	-	4.73	5.38	±	0.35	5.03	-	5.73
HGB	g/dL	6.9	±	0.4	6.5	-	7.3	13.1	±	0.6	12.5	-	13.7	16.3	±	0.8	15.5	-	17.1
	g/L	68	±	4	64	-	72	131	±	6	125	-	137	162	±	8	154	-	170
HCT	%	20.9	±	3.0	17.9	-	23.9	37.5	±	4.0	33.5	-	41.5	46.8	±	4.5	42.3	-	51.3
	L/L	0.190	±	0.030	0.160	-	0.220	0.375	±	0.040	0.340	-	0.420	0.445	±	0.045	0.400	-	0.490
MCV	fL	82.5	±	5.0	77.5	-	87.5	84.4	±	5.0	79.4	-	89.4	87.2	±	5.0	82.2	-	92.2
MCH	pg	22.5	±	3.8	18.7	-	26.3	29.7	±	3.8	25.9	-	33.5	30.2	±	3.8	26.4	-	34.0
	fmol	1.82	±	0.20	1.62	-	2.02	1.91	±	0.20	1.71	-	2.11	1.87	±	0.20	1.67	-	2.07
MCHC	g/dL	34.8	±	4.8	30.0	-	39.6	35.4	±	4.8	30.6	-	40.2	34.7	±	4.8	29.9	-	39.5
	g/L	328	±	48	280	-	376	354	±	48	306	-	402	365	±	48	317	-	413
Plt	10 ⁹ /L	82	±	25	57	-	107	221	±	45	176	-	266	461	±	65	396	-	526
MPV	fL	7.8	±	3.0	4.8	-	10.8	8.3	±	3.0	5.3	-	11.3	8.3	±	3.0	5.3	-	11.3
PCT*	%	0.02	±	0.05	-0.03	-	0.07	0.18	±	0.10	0.08	-	0.28	0.44	±	0.20	0.24	-	0.64
NEUT#	10 ⁹ /L	1.95	±	0.34	1.61	-	2.29	5.11	±	1.13	3.98	-	6.24	10.34	±	1.95	8.39	-	12.29
LYMPH#	10 ⁹ /L	1.55	±	0.47	1.08	-	2.02	1.13	±	1.00	0.13	-	2.13	2.55	±	2.36	0.19	-	4.91
MONO#	10 ⁹ /L	0.46	±	0.34	0.12	-	0.80	0.92	±	1.01	0.12	-	1.93	1.97	±	1.85	0.12	-	3.82
EO#	10 ⁹ /L	0.72	±	0.03	0.69	-	0.75	0.80	±	0.11	0.69	-	0.91	0.92	±	0.23	0.69	-	1.15
BASO#	10 ⁹ /L	0.01	±	0.01	0.00	-	0.02	0.02	±	0.02	0.00	-	0.04	0.03	±	0.03	0.00	-	0.06
NEUT%	%	57.8	±	16.6	41.2	-	74.4	63.6	±	18.4	45.2	-	82.0	66.5	±	16.0	50.5	-	82.5
LYMPH%	%	17.4	±	14.7	2.7	-	32.1	14.0	±	11.8	11.8	-	25.8	10.3	±	11.1	-0.8	-	21.4
MONO%	%	14.1	±	13.8	0.3	-	27.9	11.7	±	13.7	0.0	-	25.4	12.4	±	12.1	0.3	-	24.5
EO%	%	9.3	±	1.3	8.0	-	10.6	9.4	±	1.4	0.0	-	10.8	9.5	±	1.5	8.0	-	11.0
BASO%	%	1.5	±	0.2	1.3	-	1.7	1.4	±	0.1	0.0	-	1.5	1.4	±	0.1	1.3	-	1.5



DIAGON Ltd
Baross Str. 48-52. Budapest
Hungary-1047

