

# BC-5D

## HEMATOLOGY CONTROLS

### CONTROL

ASSAY VALUES AND EXPECTED RANGES

**LOT** BC2609B  
 2026-11-10

| Instrument                            | Parameter |                    | Low   |             | Normal |             | High  |             |
|---------------------------------------|-----------|--------------------|-------|-------------|--------|-------------|-------|-------------|
|                                       |           |                    | LOT   | BC2609BL    | LOT    | BC2609BN    | LOT   | BC2609BH    |
| BC-5800, BC-5600<br>QC Mode           | WBC       | $\times 10^9/L$    | 3.42  | $\pm 0.50$  | 7.63   | $\pm 1.00$  | 17.82 | $\pm 2.50$  |
|                                       | Neu#      | $\times 10^9/L$    | 1.53  | $\pm 0.40$  | 4.18   | $\pm 0.98$  | 11.01 | $\pm 2.15$  |
|                                       | Lym#      | $\times 10^9/L$    | 1.46  | $\pm 0.40$  | 2.56   | $\pm 0.75$  | 3.94  | $\pm 1.60$  |
|                                       | Mon#      | $\times 10^9/L$    | 0.23  | $\pm 0.23$  | 0.34   | $\pm 0.34$  | 0.93  | $\pm 0.93$  |
|                                       | Eos#      | $\times 10^9/L$    | 0.17  | $\pm 0.17$  | 0.47   | $\pm 0.47$  | 1.76  | $\pm 1.40$  |
|                                       | Bas#      | $\times 10^9/L$    | 0.03  | $\pm 0.03$  | 0.08   | $\pm 0.08$  | 0.18  | $\pm 0.18$  |
|                                       | Neu%      | %                  | 44.7  | $\pm 12.0$  | 54.9   | $\pm 12.0$  | 61.8  | $\pm 12.0$  |
|                                       | Lym%      | %                  | 42.7  | $\pm 9.0$   | 33.6   | $\pm 9.0$   | 22.1  | $\pm 9.0$   |
|                                       | Mon%      | %                  | 6.6   | $\pm 5.0$   | 4.4    | $\pm 4.4$   | 5.2   | $\pm 5.2$   |
|                                       | Eos%      | %                  | 5.0   | $\pm 5.0$   | 6.1    | $\pm 6.1$   | 9.9   | $\pm 8.0$   |
|                                       | Bas%      | %                  | 1.0   | $\pm 1.0$   | 1.0    | $\pm 1.0$   | 1.0   | $\pm 1.0$   |
|                                       | RBC       | $\times 10^{12}/L$ | 2.37  | $\pm 0.18$  | 4.09   | $\pm 0.24$  | 5.04  | $\pm 0.30$  |
|                                       | HGB       | g/L                | 66    | $\pm 4$     | 130    | $\pm 6$     | 171   | $\pm 8$     |
|                                       | HCT       | %                  | 20.2  | $\pm 2.0$   | 39.1   | $\pm 3.0$   | 52.2  | $\pm 4.0$   |
|                                       | MCV       | fL                 | 85.1  | $\pm 5.0$   | 95.6   | $\pm 5.0$   | 103.6 | $\pm 5.0$   |
|                                       | MCH       | pg                 | 27.8  | $\pm 2.5$   | 31.8   | $\pm 2.5$   | 33.9  | $\pm 2.5$   |
|                                       | MCHC      | g/L                | 327   | $\pm 30$    | 332    | $\pm 30$    | 327   | $\pm 30$    |
|                                       | RDW-CV    | %                  | 16.9  | $\pm 3.0$   | 14.1   | $\pm 3.0$   | 13.7  | $\pm 3.0$   |
|                                       | RDW-SD    | fL                 | 50.9  | $\pm 10.0$  | 48.9   | $\pm 10.0$  | 51.6  | $\pm 12.0$  |
|                                       | PLT       | $\times 10^9/L$    | 44    | $\pm 20$    | 251    | $\pm 40$    | 499   | $\pm 60$    |
|                                       | MPV       | fL                 | 8.1   | $\pm 3.0$   | 8.3    | $\pm 3.0$   | 8.9   | $\pm 3.0$   |
|                                       | PCT*      | %                  | 0.036 | $\pm 0.036$ | 0.208  | $\pm 0.100$ | 0.444 | $\pm 0.200$ |
|                                       | PDW*      | /                  | 16.9  | $\pm 3.0$   | 15.9   | $\pm 3.0$   | 16.0  | $\pm 3.0$   |
|                                       | P-LCC**   | $\times 10^9/L$    | 11    | $\pm 11$    | 47     | $\pm 25$    | 116   | $\pm 35$    |
|                                       | P-LCR**   | %                  | 24.6  | $\pm 10.0$  | 18.8   | $\pm 10.0$  | 23.3  | $\pm 10.0$  |
| BC-5390 CRP<br>BC-5310 CRP<br>QC Mode | WBC       | $\times 10^9/L$    | 3.26  | $\pm 0.50$  | 7.41   | $\pm 1.00$  | 17.06 | $\pm 2.50$  |
|                                       | Neu#      | $\times 10^9/L$    | 1.61  | $\pm 0.40$  | 4.38   | $\pm 0.98$  | 11.10 | $\pm 2.15$  |
|                                       | Lym#      | $\times 10^9/L$    | 1.30  | $\pm 0.40$  | 2.27   | $\pm 0.75$  | 3.45  | $\pm 1.60$  |
|                                       | Mon#      | $\times 10^9/L$    | 0.16  | $\pm 0.30$  | 0.26   | $\pm 0.65$  | 0.68  | $\pm 1.10$  |
|                                       | Eos#      | $\times 10^9/L$    | 0.19  | $\pm 0.20$  | 0.50   | $\pm 0.65$  | 1.83  | $\pm 1.40$  |
|                                       | Bas#      | $\times 10^9/L$    | 0.75  | $\pm 0.35$  | 1.95   | $\pm 0.80$  | 5.14  | $\pm 1.70$  |
|                                       | Neu%      | %                  | 49.6  | $\pm 12.0$  | 59.1   | $\pm 12.0$  | 65.1  | $\pm 12.0$  |
|                                       | Lym%      | %                  | 39.9  | $\pm 9.0$   | 30.7   | $\pm 9.0$   | 20.2  | $\pm 9.0$   |
|                                       | Mon%      | %                  | 4.8   | $\pm 5.0$   | 3.5    | $\pm 6.0$   | 4.0   | $\pm 6.0$   |
|                                       | Eos%      | %                  | 5.7   | $\pm 5.5$   | 6.7    | $\pm 6.5$   | 10.7  | $\pm 8.0$   |
|                                       | Bas%      | %                  | 23.0  | $\pm 10.0$  | 26.3   | $\pm 10.0$  | 30.1  | $\pm 10.0$  |
|                                       | RBC       | $\times 10^{12}/L$ | 2.31  | $\pm 0.18$  | 4.05   | $\pm 0.24$  | 4.98  | $\pm 0.30$  |
|                                       | HGB       | g/L                | 61    | $\pm 4$     | 120    | $\pm 6$     | 159   | $\pm 8$     |
|                                       | HCT       | %                  | 19.3  | $\pm 2.0$   | 38.0   | $\pm 3.0$   | 50.9  | $\pm 4.0$   |
|                                       | MCV       | fL                 | 83.7  | $\pm 5.0$   | 93.9   | $\pm 5.0$   | 102.2 | $\pm 5.0$   |
|                                       | MCH       | pg                 | 26.4  | $\pm 2.5$   | 29.6   | $\pm 2.5$   | 31.9  | $\pm 2.5$   |
|                                       | MCHC      | g/L                | 315   | $\pm 30$    | 316    | $\pm 30$    | 312   | $\pm 30$    |
|                                       | RDW-CV    | %                  | 17.2  | $\pm 3.0$   | 14.7   | $\pm 3.0$   | 14.0  | $\pm 3.0$   |
|                                       | RDW-SD    | fL                 | 50.3  | $\pm 10.0$  | 47.8   | $\pm 10.0$  | 50.1  | $\pm 12.0$  |
|                                       | PLT       | $\times 10^9/L$    | 38    | $\pm 20$    | 243    | $\pm 40$    | 492   | $\pm 60$    |
|                                       | MPV       | fL                 | 9.4   | $\pm 3.0$   | 9.1    | $\pm 3.0$   | 9.6   | $\pm 3.0$   |
|                                       | PCT*      | %                  | 0.036 | $\pm 0.050$ | 0.221  | $\pm 0.100$ | 0.472 | $\pm 0.200$ |
|                                       | PDW*      | /                  | 16.5  | $\pm 3.0$   | 16.1   | $\pm 3.0$   | 16.1  | $\pm 3.0$   |
|                                       | P-LCC**   | $\times 10^9/L$    | 9     | $\pm 15$    | 46     | $\pm 25$    | 111   | $\pm 35$    |
|                                       | P-LCR**   | %                  | 24.2  | $\pm 10.0$  | 18.9   | $\pm 10.0$  | 22.6  | $\pm 10.0$  |

# BC-5D

## HEMATOLOGY CONTROLS

### CONTROL

ASSAY VALUES AND EXPECTED RANGES

**LOT** BC2609B  
 2026-11-10

| Instrument                                                                            | Parameter |                    | Low   |             | Normal |             | High  |             |
|---------------------------------------------------------------------------------------|-----------|--------------------|-------|-------------|--------|-------------|-------|-------------|
|                                                                                       |           |                    | LOT   | BC2609BL    | LOT    | BC2609BN    | LOT   | BC2609BH    |
| BC-5300, BC-5100<br>BC-5380, BC-5180<br>QC Mode                                       | WBC       | $\times 10^9/L$    | 3.15  | $\pm 0.50$  | 7.31   | $\pm 1.00$  | 17.02 | $\pm 2.50$  |
|                                                                                       | Neu#      | $\times 10^9/L$    | 1.64  | $\pm 0.40$  | 4.42   | $\pm 0.98$  | 11.15 | $\pm 2.15$  |
|                                                                                       | Lym#      | $\times 10^9/L$    | 1.20  | $\pm 0.40$  | 2.23   | $\pm 0.75$  | 3.40  | $\pm 1.60$  |
|                                                                                       | Mon#      | $\times 10^9/L$    | 0.10  | $\pm 0.30$  | 0.09   | $\pm 0.65$  | 0.41  | $\pm 1.10$  |
|                                                                                       | Eos#      | $\times 10^9/L$    | 0.21  | $\pm 0.20$  | 0.57   | $\pm 0.65$  | 2.06  | $\pm 1.40$  |
|                                                                                       | Bas#      | $\times 10^9/L$    | 1.67  | $\pm 0.35$  | 4.66   | $\pm 0.80$  | 13.07 | $\pm 1.70$  |
|                                                                                       | Neu%      | %                  | 51.8  | $\pm 12.0$  | 60.5   | $\pm 12.0$  | 65.5  | $\pm 12.0$  |
|                                                                                       | Lym%      | %                  | 38.1  | $\pm 9.0$   | 30.5   | $\pm 9.0$   | 20.0  | $\pm 9.0$   |
|                                                                                       | Mon%      | %                  | 3.3   | $\pm 5.0$   | 1.2    | $\pm 6.0$   | 2.4   | $\pm 6.0$   |
|                                                                                       | Eos%      | %                  | 6.8   | $\pm 5.5$   | 7.8    | $\pm 6.5$   | 12.1  | $\pm 8.0$   |
|                                                                                       | Bas%      | %                  | 53.0  | $\pm 10.0$  | 63.7   | $\pm 10.0$  | 76.8  | $\pm 10.0$  |
|                                                                                       | RBC       | $\times 10^{12}/L$ | 2.35  | $\pm 0.18$  | 4.07   | $\pm 0.24$  | 5.00  | $\pm 0.30$  |
|                                                                                       | HGB       | g/L                | 62    | $\pm 4$     | 120    | $\pm 6$     | 159   | $\pm 8$     |
|                                                                                       | HCT       | %                  | 20.0  | $\pm 2.0$   | 38.9   | $\pm 3.0$   | 51.7  | $\pm 4.0$   |
|                                                                                       | MCV       | fL                 | 85.0  | $\pm 5.0$   | 95.6   | $\pm 5.0$   | 103.3 | $\pm 5.0$   |
|                                                                                       | MCH       | pg                 | 26.4  | $\pm 2.5$   | 29.5   | $\pm 2.5$   | 31.8  | $\pm 2.5$   |
|                                                                                       | MCHC      | g/L                | 310   | $\pm 30$    | 308    | $\pm 30$    | 308   | $\pm 30$    |
|                                                                                       | RDW-CV    | %                  | 16.9  | $\pm 3.0$   | 14.6   | $\pm 3.0$   | 14.2  | $\pm 3.0$   |
|                                                                                       | RDW-SD    | fL                 | 60.1  | $\pm 10.0$  | 58.1   | $\pm 10.0$  | 60.8  | $\pm 12.0$  |
|                                                                                       | PLT       | $\times 10^9/L$    | 39    | $\pm 20$    | 245    | $\pm 40$    | 494   | $\pm 60$    |
|                                                                                       | MPV       | fL                 | 8.6   | $\pm 3.0$   | 8.6    | $\pm 3.0$   | 9.0   | $\pm 3.0$   |
|                                                                                       | PCT*      | %                  | 0.034 | $\pm 0.050$ | 0.211  | $\pm 0.100$ | 0.445 | $\pm 0.200$ |
|                                                                                       | PDW*      | /                  | 16.4  | $\pm 3.0$   | 16.1   | $\pm 3.0$   | 16.1  | $\pm 3.0$   |
| BC-5000, BC-5150, Counter<br>29, BC-5120, BC-5130, BC-<br>5140, BC-5000VET<br>QC Mode | WBC       | $\times 10^9/L$    | 3.38  | $\pm 0.50$  | 7.53   | $\pm 1.00$  | 17.50 | $\pm 2.50$  |
|                                                                                       | Neu#      | $\times 10^9/L$    | 1.58  | $\pm 0.40$  | 4.35   | $\pm 0.98$  | 11.01 | $\pm 2.15$  |
|                                                                                       | Lym#      | $\times 10^9/L$    | 1.34  | $\pm 0.40$  | 2.34   | $\pm 0.75$  | 3.40  | $\pm 1.60$  |
|                                                                                       | Mon#      | $\times 10^9/L$    | 0.23  | $\pm 0.23$  | 0.24   | $\pm 0.24$  | 0.84  | $\pm 0.84$  |
|                                                                                       | Eos#      | $\times 10^9/L$    | 0.20  | $\pm 0.20$  | 0.52   | $\pm 0.52$  | 1.93  | $\pm 1.40$  |
|                                                                                       | Bas#      | $\times 10^9/L$    | 0.03  | $\pm 0.03$  | 0.08   | $\pm 0.08$  | 0.32  | $\pm 0.32$  |
|                                                                                       | Neu%      | %                  | 46.5  | $\pm 12.0$  | 57.8   | $\pm 12.0$  | 63.0  | $\pm 12.0$  |
|                                                                                       | Lym%      | %                  | 39.7  | $\pm 9.0$   | 31.1   | $\pm 9.0$   | 19.4  | $\pm 9.0$   |
|                                                                                       | Mon%      | %                  | 6.9   | $\pm 5.0$   | 3.2    | $\pm 3.2$   | 4.8   | $\pm 4.8$   |
|                                                                                       | Eos%      | %                  | 6.0   | $\pm 5.5$   | 6.9    | $\pm 6.5$   | 11.0  | $\pm 8.0$   |
|                                                                                       | Bas%      | %                  | 0.9   | $\pm 0.9$   | 1.0    | $\pm 1.0$   | 1.8   | $\pm 1.8$   |
|                                                                                       | RBC       | $\times 10^{12}/L$ | 2.35  | $\pm 0.18$  | 4.12   | $\pm 0.24$  | 5.08  | $\pm 0.30$  |
|                                                                                       | HGB       | g/L                | 62    | $\pm 4$     | 123    | $\pm 6$     | 163   | $\pm 8$     |
|                                                                                       | HCT       | %                  | 19.8  | $\pm 2.0$   | 38.7   | $\pm 3.0$   | 51.0  | $\pm 4.0$   |
|                                                                                       | MCV       | fL                 | 84.3  | $\pm 5.0$   | 93.9   | $\pm 5.0$   | 100.4 | $\pm 5.0$   |
|                                                                                       | MCH       | pg                 | 26.4  | $\pm 2.5$   | 29.9   | $\pm 2.5$   | 32.1  | $\pm 2.5$   |
|                                                                                       | MCHC      | g/L                | 313   | $\pm 30$    | 318    | $\pm 30$    | 320   | $\pm 30$    |
|                                                                                       | RDW-CV    | %                  | 20.9  | $\pm 3.0$   | 17.1   | $\pm 3.0$   | 16.5  | $\pm 3.0$   |
|                                                                                       | RDW-SD    | fL                 | 63.0  | $\pm 10.0$  | 57.7   | $\pm 10.0$  | 59.8  | $\pm 12.0$  |
|                                                                                       | PLT       | $\times 10^9/L$    | 41    | $\pm 20$    | 252    | $\pm 40$    | 507   | $\pm 60$    |
|                                                                                       | MPV       | fL                 | 10.1  | $\pm 3.0$   | 10.3   | $\pm 3.0$   | 10.9  | $\pm 3.0$   |
|                                                                                       | PCT*      | %                  | 0.041 | $\pm 0.041$ | 0.260  | $\pm 0.100$ | 0.553 | $\pm 0.200$ |
|                                                                                       | PDW*      | /                  | 15.9  | $\pm 3.0$   | 16.0   | $\pm 3.0$   | 16.3  | $\pm 3.0$   |
|                                                                                       | P-LCC**   | $\times 10^9/L$    | 12    | $\pm 12$    | 69     | $\pm 25$    | 163   | $\pm 35$    |
|                                                                                       | P-LCR**   | %                  | 29.7  | $\pm 10.0$  | 27.3   | $\pm 10.0$  | 32.2  | $\pm 10.0$  |

# BC-5D

## HEMATOLOGY CONTROLS

### CONTROL

ASSAY VALUES AND EXPECTED RANGES

**LOT** BC2609B  
2026-11-10

| Instrument                        | Parameter |                    | Low   |             | Normal |             | High  |             |
|-----------------------------------|-----------|--------------------|-------|-------------|--------|-------------|-------|-------------|
|                                   |           |                    | LOT   | BC2609BL    | LOT    | BC2609BN    | LOT   | BC2609BH    |
| BC-5390, Counter 31 AL<br>QC Mode | WBC       | $\times 10^9/L$    | 3.20  | $\pm 0.50$  | 7.40   | $\pm 1.00$  | 16.95 | $\pm 2.50$  |
|                                   | Neu#      | $\times 10^9/L$    | 1.60  | $\pm 0.40$  | 4.33   | $\pm 0.98$  | 11.02 | $\pm 2.15$  |
|                                   | Lym#      | $\times 10^9/L$    | 1.26  | $\pm 0.40$  | 2.26   | $\pm 0.75$  | 3.47  | $\pm 1.60$  |
|                                   | Mon#      | $\times 10^9/L$    | 0.16  | $\pm 0.30$  | 0.33   | $\pm 0.65$  | 0.68  | $\pm 1.10$  |
|                                   | Eos#      | $\times 10^9/L$    | 0.18  | $\pm 0.20$  | 0.48   | $\pm 0.65$  | 1.78  | $\pm 1.40$  |
|                                   | Bas#      | $\times 10^9/L$    | 0.72  | $\pm 0.35$  | 1.96   | $\pm 0.80$  | 5.09  | $\pm 1.70$  |
|                                   | Neu%      | %                  | 50.0  | $\pm 12.0$  | 58.5   | $\pm 12.0$  | 65.0  | $\pm 12.0$  |
|                                   | Lym%      | %                  | 39.5  | $\pm 9.0$   | 30.5   | $\pm 9.0$   | 20.5  | $\pm 9.0$   |
|                                   | Mon%      | %                  | 5.0   | $\pm 5.0$   | 4.5    | $\pm 6.0$   | 4.0   | $\pm 6.0$   |
|                                   | Eos%      | %                  | 5.5   | $\pm 5.5$   | 6.5    | $\pm 6.5$   | 10.5  | $\pm 8.0$   |
|                                   | Bas%      | %                  | 22.5  | $\pm 10.0$  | 26.5   | $\pm 10.0$  | 30.0  | $\pm 10.0$  |
|                                   | RBC       | $\times 10^{12}/L$ | 2.29  | $\pm 0.18$  | 4.00   | $\pm 0.24$  | 4.98  | $\pm 0.30$  |
|                                   | HGB       | g/L                | 60    | $\pm 4$     | 119    | $\pm 6$     | 158   | $\pm 8$     |
|                                   | HCT       | %                  | 19.2  | $\pm 2.0$   | 37.4   | $\pm 3.0$   | 50.5  | $\pm 4.0$   |
|                                   | MCV       | fL                 | 84.0  | $\pm 5.0$   | 93.5   | $\pm 5.0$   | 101.5 | $\pm 5.0$   |
|                                   | MCH       | pg                 | 26.2  | $\pm 2.5$   | 29.8   | $\pm 2.5$   | 31.7  | $\pm 2.5$   |
|                                   | MCHC      | g/L                | 312   | $\pm 30$    | 318    | $\pm 30$    | 313   | $\pm 30$    |
|                                   | RDW-CV    | %                  | 16.5  | $\pm 3.0$   | 14.0   | $\pm 3.0$   | 13.5  | $\pm 3.0$   |
|                                   | RDW-SD    | fL                 | 51.0  | $\pm 10.0$  | 48.0   | $\pm 10.0$  | 50.0  | $\pm 12.0$  |
|                                   | PLT       | $\times 10^9/L$    | 44    | $\pm 20$    | 245    | $\pm 40$    | 490   | $\pm 60$    |
|                                   | MPV       | fL                 | 11.2  | $\pm 3.0$   | 11.3   | $\pm 3.0$   | 11.6  | $\pm 3.0$   |
| BC-5300Vet, BC-5100Vet<br>QC Mode | WBC       | $\times 10^9/L$    | 3.15  | $\pm 0.50$  | 7.31   | $\pm 1.00$  | 17.02 | $\pm 2.50$  |
|                                   | Neu#      | $\times 10^9/L$    | 1.64  | $\pm 0.40$  | 4.42   | $\pm 0.98$  | 11.15 | $\pm 2.15$  |
|                                   | Lym#      | $\times 10^9/L$    | 1.20  | $\pm 0.40$  | 2.23   | $\pm 0.75$  | 3.40  | $\pm 1.60$  |
|                                   | Mon#      | $\times 10^9/L$    | 0.10  | $\pm 0.30$  | 0.09   | $\pm 0.65$  | 0.41  | $\pm 1.10$  |
|                                   | Eos#      | $\times 10^9/L$    | 0.21  | $\pm 0.20$  | 0.57   | $\pm 0.65$  | 2.06  | $\pm 1.40$  |
|                                   | Neu%      | %                  | 51.8  | $\pm 12.0$  | 60.5   | $\pm 12.0$  | 65.5  | $\pm 12.0$  |
|                                   | Lym%      | %                  | 38.1  | $\pm 9.0$   | 30.5   | $\pm 9.0$   | 20.0  | $\pm 9.0$   |
|                                   | Mon%      | %                  | 3.3   | $\pm 5.0$   | 1.2    | $\pm 6.0$   | 2.4   | $\pm 6.0$   |
|                                   | Eos%      | %                  | 6.8   | $\pm 5.5$   | 7.8    | $\pm 6.5$   | 12.1  | $\pm 8.0$   |
|                                   | RBC       | $\times 10^{12}/L$ | 2.35  | $\pm 0.18$  | 4.07   | $\pm 0.24$  | 5.00  | $\pm 0.30$  |
|                                   | HGB       | g/L                | 62    | $\pm 4$     | 120    | $\pm 6$     | 159   | $\pm 8$     |
|                                   | HCT       | %                  | 20.0  | $\pm 2.0$   | 38.9   | $\pm 3.0$   | 51.7  | $\pm 4.0$   |
|                                   | MCV       | fL                 | 85.0  | $\pm 5.0$   | 95.6   | $\pm 5.0$   | 103.3 | $\pm 5.0$   |
|                                   | MCH       | pg                 | 26.4  | $\pm 2.5$   | 29.5   | $\pm 2.5$   | 31.8  | $\pm 2.5$   |
|                                   | MCHC      | g/L                | 310   | $\pm 30$    | 308    | $\pm 30$    | 308   | $\pm 30$    |
|                                   | RDW-CV    | %                  | 16.9  | $\pm 3.0$   | 14.6   | $\pm 3.0$   | 14.2  | $\pm 3.0$   |
|                                   | RDW-SD    | fL                 | 60.1  | $\pm 10.0$  | 58.1   | $\pm 10.0$  | 60.8  | $\pm 12.0$  |
|                                   | PLT       | $\times 10^9/L$    | 39    | $\pm 20$    | 245    | $\pm 40$    | 494   | $\pm 60$    |
|                                   | MPV       | fL                 | 8.6   | $\pm 3.0$   | 8.6    | $\pm 3.0$   | 9.0   | $\pm 3.0$   |
|                                   | PCT*      | %                  | 0.034 | $\pm 0.050$ | 0.211  | $\pm 0.100$ | 0.445 | $\pm 0.200$ |
|                                   | PDW*      | /                  | 16.4  | $\pm 3.0$   | 16.1   | $\pm 3.0$   | 16.1  | $\pm 3.0$   |

\* For Research Use Only

\*\* These parameters are not provided on BC-5000/BC-5000 Vet analyzers  
AS904 Rev. 012