MULTI SERVICIOS LAB
INSERTOS

ASSAY SHEET

MEK-3DL HEMATOLOGY CONTROL

MEK-3DL

NIHON KOHDEN

LOT	B26YL
	2027-02-05
OPEN	14DAYS

		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305	RANGE
WBC	$10^3/\mu L$	2.1	2.1	2.09	± 0.6
RBC	$10^6/\mu L$	2.39	2.39	2.32	± 0.15
HGB	g/dL	6.0	6.1	6.14	± 0.4
HCT	%	18.5	18.5	20.4	± 2.0
MCV	fL	77.5	77.5	88.0	± 9.0
MCH	pg	25.1	25.5	26.5	± 2.4
MCHC	g/dL	32.4	33.0	30.1	± 3.0
PLT	$10^3/\mu L$	61	61	72.1	± 30
MPV	fL	5.9	8.4	8.3	± 3.0
RDW-CV *1	%	18.1	18.1	15.8	± 3.0
LY%	%	56.7	56.7	57.02	± 10.0
MO%	%	9.5	9.5	5.56	± 5.0
GR%	%	33.8	33.8	37.42	± 10.0
LY	$10^3/\mu L$	1.2	1.2	1.19	± 0.5
MO	$10^3/\mu L$	0.2	0.2	0.12	± 0.2
GR	$10^3/\mu L$	0.7	0.7	0.78	± 0.5
QC ESR *2	mm	-	-	48	± 17

*1 RDW-CV is displayed as RDW on MEK-6400, 6410, 6420, 6450.

*2 MEK-1305 only.

Manual Registration	
Lot	26Y
Level	Low
Expiration Date	2027 Y 2 M 5 D
Control	MEK-3D

MEK-1301, 1302, 1305 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1

2

3

MEK-1301, 1302, 1305 2D Barcode Entry

https://www.nihonkohden.com/spi

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN

CE 0123 IVD

0694-901684Z

Morelia (matriz) Tel: 443 690 2760
León (suc.) Tel: 477 510 2573
San Luis Potosí (suc.) Tel: 444 535 8886
www.multiservicioslab.com

MULTI SERVICIOS LAB
INSERTOS

ASSAY SHEET

MEK-3DN HEMATOLOGY CONTROL

MEK-3DN

NIHON KOHDEN

LOT	B26YN
⌛	2027-02-05
OPEN	14DAYS

		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305	RANGE
WBC	$10^3/\mu L$	7.9	7.9	7.81	± 0.8
RBC	$10^6/\mu L$	4.79	4.79	4.76	± 0.20
HGB	g/dL	13.9	14.0	14.15	± 0.5
HCT	%	41.3	41.3	46.8	± 3.0
MCV	fL	86.2	86.2	98.4	± 7.0
MCH	pg	29.0	29.2	29.7	± 2.4
MCHC	g/dL	33.7	33.9	30.2	± 3.0
PLT	$10^3/\mu L$	245	245	260.8	± 50
MPV	fL	5.4	8.1	8.2	± 2.0
RDW-CV *1	%	16.7	16.7	16.3	± 3.0
LY%	%	30.2	30.2	28.81	± 10.0
MO%	%	4.6	4.6	4.15	± 5.0
GR%	%	65.2	65.2	67.04	± 10.0
LY	$10^3/\mu L$	2.4	2.4	2.25	± 1.0
MO	$10^3/\mu L$	0.4	0.4	0.32	± 0.4
GR	$10^3/\mu L$	5.1	5.1	5.24	± 1.0
QC ESR *2	mm	-	-	9	± 5

*1 RDW-CV is displayed as RDW on MEK-6400, 6410, 6420, 6450.

*2 MEK-1305 only.

Manual Registration	
Lot	26Y
Level	Normal
Expiration Date	2027 Y 2 M 5 D
Control	MEK-3D

MEK-1301, 1302, 1305 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1

2

3

MEK-1301, 1302, 1305 2D Barcode Entry

https://www.nihonkohden.com/spi

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN

CE 0123 IVD

0694-901685Z

Morelia (matriz) Tel: 443 690 2760
León (suc.) Tel: 477 510 2573
San Luis Potosí (suc.) Tel: 444 535 8886
www.multiservicioslab.com

MULTI SERVICIOS LAB
INSERTOS

ASSAY SHEET

MEK-3DH HEMATOLOGY CONTROL

MEK-3DH

NIHON KOHDEN

LOT	B26YH
	2027-02-05
OPEN	14DAYS

		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305	RANGE
WBC	$10^3/\mu L$	21.2	21.2	20.84	± 2.2
RBC	$10^6/\mu L$	6.00	6.00	6.06	± 0.25
HGB	g/dL	19.2	19.1	19.61	± 0.7
HCT	%	55.1	55.1	63.9	± 4.0
MCV	fL	91.9	91.9	105.5	± 7.0
MCH	pg	32.0	31.8	32.4	± 2.8
MCHC	g/dL	34.8	34.7	30.7	± 3.0
PLT	$10^3/\mu L$	521	521	528.5	± 70
MPV	fL	5.0	7.9	8.0	± 3.0
RDW-CV *1	%	16.5	16.5	17.4	± 3.0
LY%	%	16.3	16.3	14.47	± 10.0
MO%	%	1.7	1.7	1.66	± 5.0
GR%	%	82.0	82.0	83.87	± 10.0
LY	$10^3/\mu L$	3.5	3.5	3.02	± 1.8
MO	$10^3/\mu L$	0.4	0.4	0.35	± 0.7
GR	$10^3/\mu L$	17.3	17.3	17.47	± 2.0

*1 RDW-CV is displayed as RDW on MEK-6400, 6410, 6420, 6450.

Manual Registration	
Lot	26Y
Level	High
Expiration Date	2027 Y 2 M 5 D
Control	MEK-3D

MEK-1301, 1302, 1305 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1

2

3

MEK-1301, 1302, 1305 2D Barcode Entry

https://www.nihonkohden.com/spi

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN

CE 0123 IVD

0694-901686Z

Morelia (matriz) Tel: 443 690 2760
León (suc.) Tel: 477 510 2573
San Luis Potosí (suc.) Tel: 444 535 8886
www.multiservicioslab.com